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1. Atrofia de iris

Fluorescein angiographic studies by *Strobel and Izák* have shown no leakage of the iris vessels at the enclavation sites. In addition no iris atrophy has been reported in the area of the fixation. Since its invention 25 years ago, the fixation principle of the ARTISAN® Lens has remained unchanged.

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

3. Daño endotelial

a) Respetar la ACD mínima requerida
3.0 mm desde epitelio > Artisan

¹Tahzib NG, Nuijts RM, Wu WY, Budo CJ Long-term Study of Artisan Phakic Intraocular Lens Implantation for the Correction of Moderate to High Myopia Ten-Year Follow-up Results. Ophthalmology 2007; 114(6):1133-42

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



¹ Posted in Feb 2016 · Cataract and Refractive · Eurotimes Stories · <http://www.eurotimes.org/node/2266>

LOOMING EPIDEMIC by Roibeard O'hEineachain¹

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1. Retrospective Study of the ARTISAN Aphakia IOL with 10 year follow-up

Nayyirih G. Tahzib, Ophthalmology
2007; 114(6):1133-42

2. Paired-eye comparison of corneal endothelial cell counts after unilateral iris-claw phakic intraocular lens.

Merce Morral, MD, PhD, José L. Güell, MD, PhD, Mostafa A. El Hussein, MD, Daniel Elies, MD, Oscar Gris, MD, PhD, Felicidad Manero, MD Journal of Cataract & Refractive Surgery January 2016

• Long-term endothelial cell loss with the iris-claw intraocular Phakic lenses (Artisan®)

Virgilio Galvis^{1,2,3} & John F. Villamil^{1,4} & María Fernanda Acuña^{2,3} & Paul A. Camacho^{2,3} & Jesús Merayo-Llaves^{5,6} & Alejandro Tello^{1,2,3} & Sandra Lizeth Zambrano^{2,4} & Juan José Rey³ & Juan Vicente Espinoza^{2,4} & Angélica María Prada¹

– Received: 25 June 2019 /Revised: 19 August 2019 /Accepted: 7 October 2019

Long-term Study of Artisan Phakic Intraocular Lens Implantation for the Correction of Moderate to High Myopia

Ten-Year Follow-up Results

Nayyirih G. Tahzib, MD,¹ Rudy M. Najib, MD, PhD,² Wensly Y. Wu, MD,¹ Camille J. Budo, MD¹

Purpose: To determine the long-term performance of the Artisan phakic intraocular lens (PIOL) for the correction of myopia.

Design: Long-term (10 years) retrospective follow-up study.

Participants: Eighty-nine eyes of 49 patients who underwent Artisan PIOL implantation for the correction of myopia.

Methods: Comparisons were made between preoperative clinical data and the clinical data at 1, 6, and 10 years after surgery.

Main Outcome Measures: Refractive stability, refractive predictability, safety, efficacy, best-corrected visual acuity (BCVA), uncorrected visual acuity (UCVA), intraocular pressure, intraoperative problems, corneal endothelial cell density, corneal endothelial cell loss, and glare levels were evaluated.

Results: The mean spherical equivalent (SE) after 10 years was -0.70 ± 1.00 diopters (D; range, -4.00 to 2.00 D), with no significant change in mean SE between 1, 6, and 10 years. At 10 years, 68.8% of all eyes were within 1.0 D of the intended correction. At 10 years, 31.2% ($n = 24$) gained 1 or more Snellen lines of BCVA and 2.6% ($n = 2$) lost more than 2 Snellen lines of BCVA; 93.3% reached a BCVA of 20/40 or better, and 82.0% reached a UCVA of 20/40 or better. The mean intraocular pressure remained stable and was 15.5 ± 3.5 mmHg (range, 7–25 mmHg) at 10 years. The mean endothelial cell loss was $-8.86 \pm 16.01\%$ (range, -51.69% to 34.43%) at 10 years.

Conclusions: Long-term results demonstrate that the implantation of an Artisan PIOL for the correction of moderate to high myopia is a stable, predictable, and safe method when strict inclusion criteria for surgery are applied. There was no significant loss of corneal endothelial cells and no reports of long-term glare. *Ophthalmology* 2007;114:1133–1142 © 2007 by the American Academy of Ophthalmology.

Gratias Actio for Clinical and Experimental Ophthalmology
https://doi.org/10.1007/s00417-019-04069-9

REFRACTIVE SURGERY

Long-term endothelial cell loss with the iris-claw intraocular phakic lenses (Artisan®)

Virgilio Galvis^{1,2,3} · John F. Villamil^{1,4} · María Fernanda Acuña^{2,3} · Paul A. Camacho^{2,3} · Jesús Merayo-Llaves^{5,6} · Alejandro Tello^{1,2,3} · Sandra Lizeth Zambrano^{2,4} · Juan José Rey³ · Juan Vicente Espinoza^{2,4} · Angélica María Prada¹

Received: 25 June 2019 / Revised: 19 August 2019 / Accepted: 7 October 2019
© Springer-Verlag GmbH Germany, part of Springer Nature 2019

Abstract

Purpose: To evaluate the endothelial cell loss in patients with iris-claw phakic lenses (Artisan®) in a long-term follow-up.

Methods: We analyzed the medical records of patients who had undergone iris-claw phakic lens implantation and who had at least 5 years of follow-up.

Results: We included 67 eyes with myopic errors (follow-up 9.6 ± 3.0 years) and 10 eyes with mixed astigmatism or hyperopic errors (follow-up 8.8 ± 2.5 years). The mean total endothelial density loss at the last follow-up visit was $18.5\% \pm 17.0\%$ and $10.5\% \pm 12.3\%$, respectively. 29.9% of the eyes in the myopic group and 20% in the hyperopic group lost more than 25% of the preoperative endothelial cell density. During the postoperative follow-up period, 60.8% of the eyes in the myopic group and 40% of the eyes in the hyperopic group lost a higher percentage than the expected physiological loss. Two eyes in the myopic group (3.0%) had a final cell density of less than 1200 cells/mm². None of the variables studied had a statistically significant association with the percentage of annual endothelial loss in the postoperative period. Three phakic lenses were explanted: two by cataract and one by contact and severe decrease of the endothelial density (862 cells/mm²).

Conclusions: There is a significant endothelial cell loss in a low percentage of the eyes with Artisan® lenses in the long term, and it can decrease to critical levels. Periodic endothelial density evaluations are required for these patients. The selection criteria of surgical candidates could be reevaluated.

Keywords: Phakic intraocular lenses · Corneal endothelial cell loss · Corneal edema · Refractive surgery · Corneal endothelium · Cornea

Introduction

LASIK surgery has limitations for the correction of very high refractive errors (approximately above 8.00 Diopters (D)) for myopia and 5.00 D for hyperopia), since considerable undesirable effects of the photorefractive procedure may occur, such as loss of vision, decreased contrast sensitivity, and night vision symptoms, as well as long-term regression of the refractive correction. In addition, the excessive corneal tissue ablation necessary for these high corrections is a risk factor for iatrogenic corneal ectasia [1–3]. As an alternative in these cases, various intraocular lenses have been designed to be implanted in these patients without removing the lens, known as phakic intraocular lenses (PIOLs) [4–7]. Within this group are the angle-supported lenses that were initially implanted in the mid-1950s and withdrawn a few years later due to serious complications. New

✉ Alejandro Tello
alejandro@grad.ums

¹ Centro Oftalmológico Virgilio Galvis, FOSCAL Internacional, Calle 108 30-95, Condesado 301, Tere C, Floridablanca, Colombia

² Ophthalmology Department, Fundación Oftalmológica de Santander FOSCAL, Floridablanca, Colombia

³ Faculty of Health Sciences, Universidad Autónoma de Bucaramanga, Bucaramanga, Colombia

⁴ Faculty of Health Sciences, Universidad Industrial de Santander, Bucaramanga, Colombia

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⁶ Ophthalmology Department, Universidad de Oviedo, Oviedo, Spain

ARTISAN & ARTIFLEX Phakic IOL [Training Manual]

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

Use of a security suture during retropupillary implantation of an iris-claw IOL, Joshua C. Teichman, MD

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The screenshot shows a YouTube channel page for 'Ophtecbv'. The channel name is 'Artisan Aphakia' with 22 videos and 1,391 views. The page displays two columns of video thumbnails and titles. The left column shows a playlist of 10 videos, and the right column shows a list of 10 individual videos. Each video title includes the name of the procedure and the surgeon's name.

Video Number	Video Title	Channel
1	ARTISAN Aphakia IOL in Pediatric surgery by Dr. J.L. Güell	Ophtecbv
2	Retro-pupillary ARTISAN in patient with Axenfeld-Rieger and Marfan's by Dr. A. John Kanelopoulos	Ophtecbv
3	ARTISAN Aphakia implantation using the enclavation forceps Dr C. Peckar	Ophtecbv
4	ARTISAN Aphakia "The ideal IOL for secondary implantation" by Dr. J.L. Güell	Ophtecbv
5	ARTISAN Aphakia Retropupillar "The Dark side of the moon" implantation by Dr. A. Mohr	Ophtecbv
6	Artisan Aphakia IOL implantation by dr. Adrián Hernández Martínez	Ophtecbv
7	Artisan Aphakia implantation by Sergio Bonafonte, Barcelona	Ophtecbv
8	Artisan Aphakia enclavation using VacuFix by Dr. V.I. Apostolov	Ophtecbv
9	Femto DSAEK and ARTISAN IOL exchange by Dr. Kanelopoulos	Ophtecbv
10	Artisan Aphakia retro-pupillary by Dr. M Marcos Robles	Ophtecbv
1	Enclavation of The Artisan PIOL by Budo	Ophtecbv
2	ARTISAN PIOL implantation technique by Dr. C. Budo	Ophtecbv
3	ARTIFLEX PIOL implantation technique by Dr. C.Budo	Ophtecbv
4	Needle enclavation of the ARTIFLEX PIOL	Ophtecbv
5	ARTIFLEX implantation. Lateral. By Dr J.L. Güell	Ophtecbv
6	ARTIFLEX enclavation using the VacuFix by Dr. J.L. Güell	Ophtecbv
7	ARTIFLEX implantation in depth Part 1, Enclavation pointing by Chan Young Im	Ophtecbv
8	ArtiLens radio commercial	Ophtecbv
9	Vacufix vs Enclavation needle: Comparative study By: Dr. O. BARRADA, Egypt	Ophtecbv
10	ARTIFLEX implantation by Dr J.L. Güell	Ophtecbv

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<https://www.ophtec.com/>

<https://www.youtube.com/user/Ophtecbv>